
Stereo microscope Stemi 305 cam with integrated camera

**Operating manual
Bedienungsanleitung
Mode d'emploi
Instrucciones de manejo
Istruzioni per l'uso
Manual de instruções
操作手册
取扱説明書**

English

Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

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Issued by: Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena, Germany
microscopy@zeiss.com
www.zeiss.com/microscopy

Carl Zeiss Microscopy GmbH
Königsallee 9-21
37081 Göttingen, Germany

Number of this manual: 435063-9044-008
Date of issue: Version 5, 03/20/2015

CONTENTS

	Page
1 INTRODUCTION	2
1.1 General information	2
1.2 Notes on Device Safety	3
1.3 Warranty Notes	7
2 DESCRIPTION	8
2.1 Indication for use	8
2.2 Microscope system	8
2.3 Interfaces on the microscope body Stemi 305 cam	9
2.4 System Overview	10
2.5 Technical Data	12
3 INSTALLATION	16
3.1 General information	16
3.2 Connecting stereo microscope to the power supply	16
3.3 Switching stereo microscope ON or OFF	16
4 OPERATION OF THE INTEGRATED CAMERA WITH WI-FI.....	17
4.1 Commissioning the camera	17
4.1.1 Mode "WLAN Access-Point" (factory setting)	18
4.1.2 Mode "Connection to existing WLAN"	19
4.2 Setting the WLAN mode (Wi-Fi) of the integrated camera	19
4.3 Reset feature	21
5 CARE, MAINTENANCE AND SERVICE	22
5.1 Care	22
5.2 Maintenance	23
5.2.1 Wearing parts	23
5.3 Service	23
5.4 Product disposal	23
6 APPENDIX	24
6.1 List of abbreviations	24
6.2 Troubleshooting	25
6.3 List of illustrations	26
6.4 Index	27

1 INTRODUCTION

1.1 General information

The stereo microscope Stemi 305 cam has been designed, produced and tested in compliance with the standards DIN EN 61010-1 (IEC 61010-1) and IEC 61010-2-101 "Safety Requirements for Electrical Measuring, Control and Laboratory Equipment".

The power unit used corresponds to DIN EN 60950-1 (IEC 60950-1) "Information technology equipment - Safety".

This device complies with the requirements of the EC Directive 98/79/EC Annex 1 regarding ivd products, EC RoHS Directive 2011/65/EG (RoHS) and EC Directive 1999/5/EG (R&TTE, Radio Equipment & Telecommunications Terminal Equipment) and bears the **CE** - mark.

This operating manual contains information and warnings to be observed by the operator.

The devices must be disposed of in compliance with the WEEE Directive 2012/19/EC.

The warning and information symbols explained below are used in these operating instructions:

**CAUTION**

This symbol identifies a risk for the user.

**CAUTION**

*LED risk group 2 acc. to DIN EN 62471:2009 Optical radiation is emitted.
Do not look into the beam. This may damage your eyes.*

**CAUTION: High-energy UV radiation!**

Risk of damage to eyes and skin!

**CAUTION**

Hot surface!

**CAUTION**

Remove the mains plug before any intervention in the instrument!

**ATTENTION**

This symbol identifies a risk for the device or system.

**NOTICE**

This symbol identifies a note which must be given special attention.

1.2 Notes on Device Safety

If it becomes apparent that the safety measures are no longer effective, the device must be taken out of service and secured against being used again unintentionally. Please contact ZEISS customer service or the Carl Zeiss Microscopy Service to have the device repaired.

The Stemi 305 Information regarding device safety and installation conditions 435063-6044-008 and the Stemi 305 Concise operating instructions 435063-8044-008 as well as the operating instructions regarding the light sources also must be observed.

Do not operate the devices included in the scope of supplies in a potentially explosive atmosphere, in the presence of volatile anaesthetics or flammable solvents such as alcohol, gasoline or similar substances.

The devices may only be operated by instructed persons who have been informed about the possible risks in conjunction with the use of microscopes and the fields of application in question. The microscope is a high-precision device that may be impaired in its performance or even destroyed if handled improperly.

Before commissioning the device, make sure that the existing power supply is suitable for the device.

Make sure that the mains plug is removed from the socket before the device is opened.

The devices are not equipped with special equipment protecting them from corrosive, potentially infectious, toxic and radioactive or other samples that may be hazardous to health. If you handle such samples, be sure to observe all legal requirements, in particular the relevant national accident prevention regulations.

The LED reflected light illuminators and the LED transmitted light illuminators have been classified as belonging to LED risk group 2 acc. to DIN EN 62471:2009. Avoid looking directly into the LED light.

If the device is operated with an external cold light source (high-energy light), never look directly into the fibre optic output of the cold light source. Otherwise, there is a risk of dazzling and blindness.

The Stemi 305 cam body with Wi-Fi may only be used if it is covered by a radio license for your region. For further information on the radio licences, see Annex to the operating manual. If necessary, contact your regional ZEISS representative.

The visualized images of the Stemi 305 cam body with Wi-Fi may only be used for professional training and research. These images are not meant to be used for direct generation of diagnostic results.

Never cover the open light guide socket or the fibre optic output. There is a fire hazard!

In all cases avoid covering the open light guide socket or the fibre optic output with your hand or other parts of your body. There is a risk of burns!

Safe disconnection from the power supply is ensured exclusively by removing the mains plug. The switch on the microscope only switches into standby mode.

Do not replace detachable mains cables by inadequately designed mains cables. Only the specified power supply lines may be used.

The microscope and the Controller K LED each feature a desktop power supply which permits changing the voltage setting on the instrument.

The desktop power unit meets the requirements of protection class II (with protective insulation). If its casing is damaged, put the power unit out of operation. The microscope may only be operated with the desktop power unit supplied with the instrument.

Changing the lamps on the cold light source is only admissible in accordance with the manufacturer's operating manual. If the appropriate instructions are not complied with, there is a risk of burns and explosion when changing the lamp.

Defective devices should not be disposed of with household waste; dispose of them in compliance with the statutory provisions.

Samples must also be disposed of appropriately in accordance with the applicable statutory provisions and internal work instructions.

The stereo microscope Stemi 305 cam incl. its original accessories may only be used for the applications described in this operating manual. The manufacturer cannot assume any liability for other applications, including those of individual modules or single components.

Changes to and maintenance on this device and on devices which are operated together with the microscope may only be performed by our service or by authorized persons. The manufacturer of the devices does not accept any liability for damage due to unauthorized interventions in the device. Moreover, such interventions invalidate all guarantee / warranty claims.

Dirt and dust may impair the performance of the device. The device must therefore be protected from such influences to the greatest possible extent and covered with the dust cover when not in use. Before covering the device, always make sure that it is definitely switched OFF. Considerable temperature fluctuations, direct exposure to sunlight and vibrations should be avoided.

Blocking or covering ventilation louvers may result in heat accumulation which may damage the device and in extreme cases even cause fire. Always keep ventilation louvers free and do not insert or allow any objects to fall into them. All electrical components and parts must be located at a distance of min. 15 cm from flammable objects and walls.

If it is intended to operate the Stemi 305 cam with an external fibre optical cold light source, please refer to the appropriate operating instructions incl. the safety provisions contained in it before commissioning the light source.

The desktop power unit must not get into contact with moisture.

The Stemi 305 cam body with Wi-Fi may only be used if it is covered by a radio license for your region. Contact your ZEISS representative.

To protect the WLAN against access by unauthorized persons while operating the Stemi 305 cam, it is recommended to change the standard password via the WLAN configuration website. When selecting the new WLAN password, make sure that the password contains digits and letters and is as long as possible.

Changes to the Stemi 305 cam body involving Wi-Fi (e. g. aerials) are not admissible, as this would invalidate the existing radio licenses.

Should the WLAN password for the Stemi 305 cam no longer be known, the device can be reset to the factory settings via the Reset feature (incl. WLAN-SSID and WLAN password, see data on the back of the device).

For transport over an extended distance, the device must be dismantled partially and transported in its original packaging.

To transport the device over an extended distance, use the original packaging or the transport case Stemi 305/508.



Fig. 1 Warning and information labels on Stemi 305 cam

- 1 LED aperture, transmitted-light illuminators
- 2 LED aperture, spotlight K LED
- 3 LED aperture, ring light K LED, segmentable
- 4 LED aperture, double spotlight K LED

Fig. 2 Warning labels and apertures for LED radiation

1.3 Warranty Notes

The manufacturer guarantees that the device is free from material or manufacturing defects when delivered. Any defects which may have occurred must be reported to us immediately and steps be taken to minimize damage. If notified of such a defect, the device manufacturer is obligated to rectify it at his discretion, either by repairing the instrument or by delivering an intact replacement. No guarantee is provided for defects caused by natural wear (wearing parts in particular) and improper use.

The instrument manufacturer shall not be liable for damage caused by faulty operation, negligence or any other tampering with the device, particularly the removal or replacement of instrument components, or the use of accessories from other manufacturers. Such actions will render all warranty claims invalid.

No maintenance or repair work may be performed on the microscopes which exceed the activities specified in this operating manual. Repair may only be performed by ZEISS Service or by persons specifically authorized by it. Should any malfunctions occur on the device, please first contact the Carl Zeiss Microscopy Service or, abroad, the ZEISS representative in your area.

2 DESCRIPTION


2.1 Indication for use

The stereo microscope Stemi 305 cam is a microscope for magnifying spatial observation of small objects, visualization of the objects and wireless transmission of the images via WLAN. It has been designed and built for training purposes in schools, universities and natural science facilities. It is used, moreover, in biological and medical laboratories, in industrial manufacture and quality assurance.

Stemi 305 cam is envisaged for applications in biology and in medicine for the analysis of blood and / or tissue samples from the human body. Applications in the field of diagnostic medicine are explicitly excluded, except for the field of medical research.

The visualized images of Stemi 305 cam must be used for training, laboratory and research purposes. These images are not meant to be used for direct generation of diagnostic results.

2.2 Microscope system

- | | | | |
|---|---|----|--|
| 1 | 1 | 8 | Rotary/push-button for switching ON and OFF the transmitted-light illuminator and for adjusting the illumination intensity |
| 2 | 2 | 9 | Reflected-light illuminator (spotlight K LED) |
| 3 | 3 | 10 | Lever for adjusting the transmitted-light illuminator – bright field or dark field on the stand K EDU |
| 4 | 4 | 11 | Insert plate for placing the specimen |
| 5 | 5 | | |
| 6 | 6 | | |
| 7 | 7 | | |

Fig. 3 Microscope system Stemi 305 cam

The illustrated microscope equipment shows an example and may differ from the ones actually existing!

2.3 Interfaces on the microscope body Stemi 305 cam

- 1** RJ-12 socket:
For the power supply of the integrated camera and the integrated vertical illumination
- 2** Status LED:
Indicates the camera's operating state:
 - flashes while initialization is in progress (the camera starts 3 or 4 seconds after switching-ON),
 - is permanently lit when the camera is ready to operate.
- 3** Reset opening:
To reset the integrated camera




Fig. 4 Interfaces on the Stemi 305 cam body


2.4 System Overview

English


Stemi 305 MAT microscope set
435063-9030-100


Stemi 305 EDU microscope set
435063-9010-100


Stemi 305 LAB microscope set
435063-9020-100


Stand K MAT
435425-9020-000
(incl. plastic plate B/W
and dust cover)


Stand K EDU
435425-9000-000
(incl. glass plate
and dust cover)


Stand K LAB
435425-9010-000
(incl. glass plate
and dust cover)

Stand K
435424-9000-000
(incl. plastic plate B/W
and dust cover)

Polarizer transillumination K/M
435526-9020-000

Double spot illuminator K LED
435525-9020-000

Spot illuminator K LED
435525-9010-000

Polarizer spot K LED
435526-9010-000

as accessory:
Hand rest stand K LAB
435425-9010-010

as a substitute:
Glass plate, D=84x5 mm
435425-9310-000

as a substitute:
Plastic plate B/W, D=84x5 mm
435425-9320-000


as a substitute:
Dust cover
415500-1800-000

2.5 Technical Data

Stemi 305 cam in stand K EDU

Dimensions


Weight

Stemi 305 cam in stand K EDU..... 4.6 kg

Stemi 305 cam in stand K LAB

Dimensions

Weight

Stemi 305 cam in stand K LAB 6.2 kg

Ambient conditions

Storage (in packaging)

Admissible ambient temperature +10 °C to +40 °C
 Admissible relative humidity max. 75 % at +35 °C (non-condensing)

Transport (in packaging)

Admissible ambient temperature -40 °C to +70 °C

Operation

Admissible ambient temperature +10 °C to +40 °C
 Admissible relative humidity max. 75 %
 Atmospheric pressure 800 hPa to 1060 hPa
 Pollution degree..... 2
 Field of application..... closed rooms
 Height of field of application max. 2000 m

Operating data – desktop power units, microscope and controller K LED

Safety Class II
 Type of enclosure IP 20
 Electrical safety acc. to DIN EN 61010-1 (IEC 61010-1)
 taking CSA and UL regulations into account
 Pollution degree..... 2
 Overvoltage category 2
 Line voltage 100 V to 240 V ±10 %
 Due to the fact that the instrument features a multi-voltage power supply unit, change-over of the device voltage not required!
 Line frequency 50 Hz – 60 Hz
 Power consumption: Desktop power unit with connected microscope max. 40 VA
 Output desktop power unit for microscope
 and for Controller K LED 12 V DC, max. 2 A
 Input desktop power unit for microscope
 and for Controller K LED 100 V to 240 V, 50 Hz – 60 Hz, max. 0.55 A

Operating data – Stemi 305 cam (Integrated HD IP Wi-Fi camera 1.2 MP with integrated camera adapter 0.5x)

Sensor	1/3" colour, CMOS
Resolution	1600 (H) x 1200 (V) = 2.0 megapixels
Pixel size	2.8 µm x 2.8 µm
Active sensor area	4.73 mm x 3.52 mm
All images are scaled down to 1.2 MP for Wi-Fi streaming.	
Single image (downscaled)	1280 pixels x 960 pixels (JPEG)
IP livestream (downscaled)	1280 pixels x 960 pixels @ 15 fps (MJPEG)
Dynamic range	71 dB
Signal-to-noise ratio	42.3 dB
Interface	Wi-Fi (IEEE 802.11b/g/n, 2.4 GHz, 6 dBm)
WLAN range	5 m – 15 m, depending on viewing direction and number of WLAN/Bluetooth equipment items within the range
White balance	AWB and various white balance modes
Frame rate	15 fps
Wi-Fi-Standard	IEEE802.11b/g/n
Wi-Fi output power	6 dBm
Access format	Infra mode
SSID	configurable by user
Bandwidth	2.4 GHz
Hardware:	
Stemi 305 cam	Version 1.0 or above
Tablet PC	Apple iPad Version 4 or above
Desktop power unit	Model GFP241DA-1220B-1 or more recent
Software	
WI-FI camera	Firmware version: E20141007 or higher
iPad apps	Labscope and Matscope Version 1.1 or higher
Controllable via	two different WLAN modes

1. Mode "WLAN Access Point" (factory setting)
Without additional WLAN equipment, up to 6 iPads can be connected to **one** microscope with Stemi 305 cam body – in this mode, the microscope provides the WLAN (standard mode).
2. Mode "Connection to existing WLAN"
One or several microscopes with Stemi 305 cam body can be integrated in an existing WLAN – in this mode, **one** iPad can "see" all microscopes simultaneously. Change-over between the various microscope live images is easy and comfortable.

The Stemi 305 cam body has been certified for sale in the following countries: EU, USA, Japan, China. Certification for further countries is under way. In case of doubt, contact the appropriate ZEISS representative.

Optical risk group classification acc. to DIN EN 62471:2009

Overall device	LED risk group 2 acc. to DIN EN 62471:2009
Integrated vertical illumination	LED risk group 2 acc. to DIN EN 62471:2009
Spotlight K LED (reflected illumination)	LED risk group 2 acc. to DIN EN 62471:2009
Double spotlight K LED (reflected illumination).....	LED risk group 2 acc. to DIN EN 62471:2009
Transmitted light unit in stand K LAB	LED risk group 2 acc. to DIN EN 62471:2009
Transmitted light unit in stand K EDU	LED risk group 2 acc. to DIN EN 62471:2009
LED reflected illumination, peak at 460 nm,	LED risk group 2 acc. to DIN EN 62471:2009
LED transmitted illumination, peak at 465 nm,	LED risk group 2 acc. to DIN EN 62471:2009

3 INSTALLATION

3.1 General information

-  Before installing and commissioning the device, make sure to read the **Notes on Device Safety** carefully (see section 1.2, Page 3).

The Stemi 305 cam with the necessary tools and optional equipment is delivered in several standard packages.

- Take all units out of the packaging and check them for completeness according to the delivery note.

-  For installation and / or transport of the microscope or the stand, exclusively use the handle (Fig. 3/4) provided to this effect.

- Remove transport securing devices (adhesive tapes etc.).
- Keep original packaging for a possibly extended period of storage or return of the device to the manufacturer or dispose of it appropriately.


Fig. 5 Connecting the stereo microscope

3.2 Connecting stereo microscope to the power supply

- Insert the mains cable (Fig. 5/2) in the mains supply socket (Fig. 5/1) of the stand.
- Connect the mains cable (Fig. 5/2) to a mains supply socket outlet.

3.3 Switching stereo microscope ON or OFF

- Use the power switch (Fig. 5/3) to switch the stereo microscope ON or OFF.

-  Utilization and configuration of the Stemi 305 cam body (with integrated camera) has been described in section 4, as of page 17.

4 OPERATION OF THE INTEGRATED CAMERA WITH WI-FI

The stereo microscope is placed on a level, stable base, installed acc. to the operating manual Stemi 305 (435063-7044-010), Chapter 3 and connected to the power supply system.

The basic setting of the microscope and the illuminators has been effected previously acc. to Chapter 4 of the operating manual Stemi 305.

An Apple iPad (to be acquired separately from electronics retailers) is available and has been switched ON.

The ZEISS App Labscope or Matscope have been downloaded from the Apple Appstore and installed on the iPad.

System requirements:

- Apple iPad 4 or more recent
- Labscope / Matscope App Version 1.1 or higher

4.1 Commissioning the camera

- Switch the stereo microscope ON acc. to section 3.3 using the power switch (Fig. 5/3).

After a few seconds, the blue status LED in the microscope body (Fig. 6/1) starts flashing, which indicates that initialization of the integrated Wi-Fi camera is in process.

- Wait until the blue status LED is lit permanently.

On initial commissioning, the camera is in the mode "WLAN Access-Point" (factory setting), i.e. after initialization of the camera, the microscope provides a WLAN.

To use the installed camera, an iPad must be connected to this WLAN.

- To this effect, select in the iPad settings the SSID of the Stemi 305 cam network and enter the appropriate password.
- Subsequently, the apps Labscope or Matscope can be started which recognize the camera automatically. There, the live image appears as thumbnail. To continue using the Labscope / Matscope, follow the information screen of the apps.


Fig. 6 Status LED in microscope body


Fig. 7 Adhesive label on the back of the device

- The adhesive label on the back of the microscope body (Fig. 7) shows the WLAN identifier (SSID) set by the factory, and the WLAN password of the device. SSID and WLAN password can be changed subsequently via a WLAN configuration website.

SSID: **Stemi-xyyyzz**
WLAN password: **ZEISS1846 [...]**

At choice, the mode "Connection to existing WLAN" can be used for camera control.

The Stemi 305 cam body with Wi-Fi may only be used if it is covered by a radio license for your region. Contact your ZEISS representative.

The visualized images of the Stemi 305 cam body with Wi-Fi may only be used for professional training and research. These images are not intended to be used for direct generation of diagnostic results.

4.1.1 Mode "WLAN Access-Point" (factory setting)

In this mode, each Stemi 305 cam provides its own WLAN.

In the WLAN settings of the iPad, all Stemi 305 WLAN networks are recognized which are within reach.

After establishing the connection to one of these WLAN networks, the connected Stemi 305 cam appears as an icon in the Apps Labscope / Matscope and can be controlled via the iPad. Change-over to another Stemi 305 cam is only possible if the user changes over to the WLAN of the required Stemi 305 cam in the WLAN settings of the iPad.

Depending on the WLAN and on the image contents, between 3 and max. 6 iPads can be connected to the WLAN of a Stemi 305 cam, thus being used for viewing live images, for camera control or snapshots.

Due to the delayed live stream and the line loss thus caused, it is not recommended to use more than 6 live streams.

To operate several Stemi 305 cams comfortably (without changing the WLAN network), each Stemi 305 cam disposes of another operating mode, i.e. the mode "Connection to existing WLAN" which can be set via a WLAN configuration website.

4.1.2 Mode "Connection to existing WLAN"

In this mode, after entry of the password, the Stemi 305 cam connects to an existing WLAN (can be selected by the user).

It can be provided either by another Stemi 305 cam or by a separate WLAN router.

Upon connecting several Stemi 305 cams and one or several iPads to the same WLAN, all connected microscopes are displayed simultaneously in the App Labscope / Matscope. Thus, the user can change between the microscopes with ease.

When using the WLAN of one of the Stemi 305 cams, approx. 3 microscopes and 3 iPads can be interconnected without this delaying the live stream.

Via an external WLAN router of the N300 class, 10 microscopes and 10 iPads can be used without problems. If additional (connected via cables) WLAN Access Points are used, the WLAN can be scaled up even further. Thus, scenarios with multiple microscopes / iPads can be implemented in the same tutorial room.

4.2 Setting the WLAN mode (Wi-Fi) of the integrated camera

The WLAN mode of the Stemi 305 cam can be changed as follows by means of the configuration website (**Configure Microscope**):

- Start the App Labscope on the iPad (see Fig. 8).


Fig. 8 Configuration page of Stemi 305 cam

- Touch the button **Advanced camera settings** and, in the following pop-up window, the link **WLAN configuration web-interface**.

The Safari browser of the iPad opens (see Fig. 9).


Fig. 9 Safari browser of the iPad – selection of WLAN mode

- Select the required mode **As an access point** (Mode "WLAN Access Point") or **Connect to existing WLAN** (Mode "Connection to existing WLAN") and acknowledge by pressing **Next**.
- Close Safari Browser and configuration page.

If there is no iPad, and provided the Stemi 305 cam is in mode "WLAN Access Point", the following approach is possible as an alternative:

- Connect PC to the WLAN of the Stemi 305 cam and open an available browser on the PC.
- Enter **http://192.168.1.26** into the address line.
- Enter user name **admin** and password **ZEISS1864**.

The configuration page **Configure Microscope** in the browser opens.

- Change over the WLAN setting as described above to the mode "Connection to existing WLAN".

If the Stemi 305 cam is in mode "Connection to existing WLAN", the mode cannot be changed directly. To this effect, the Stemi 305 cam must be reset to the factory settings via the Reset feature.

4.3 Reset feature

The Wi-Fi camera may require resetting if the integrated camera

- is in mode "Connection to existing WLAN", but there is no WLAN available,
- is in mode "WLAN Access Point", but the user forgot the WLAN password,
- behaves atypically.

After the reset to the factory settings, the microscope camera is in mode "WLAN Access Point".

For the WLAN identifier (SSID and password), please refer to the label on the back of the device (see also section 4.1, page 17).


Fig. 10 Reset of the integrated camera

To reset the Wi-Fi camera, proceed as follows:

- Switch Stemi 305 cam OFF and ON again.
- Wait for camera initialization to be completed.
- Keep the Reset button (Fig. 10/1) pressed for min. 10 seconds using an unbent paper clip.
- Wait for the Reset procedure to be completed.
This takes approx. 80 to 100 seconds.

During the Reset procedure, the status LED of the camera first goes out, then flashes twice. Camera initialization being complete, the status LED (Fig. 10/2) is lit permanently.

5 CARE, MAINTENANCE AND SERVICE

5.1 Care

Care of the devices is restricted to the following operations:

The devices are not equipped with special equipment protecting them from corrosive, potentially infectious, toxic and radioactive or other samples that may be hazardous to health. If you handle such samples, be sure to observe all legal requirements, in particular the relevant national accident prevention regulations.

- Remove contamination in accordance with the national accident prevention regulations.
- Disconnect the devices from the mains power after use. Protect the devices from dust and moisture by using a suitable cover (protective device cover).
- Never expose the devices to inadmissible climate conditions (high humidity and temperature).

Disconnect the devices from the mains power before cleaning. Make sure, that no cleaning fluids penetrate into the device.

Tenacious impurities on glass surfaces such as fingerprints and traces of grease, are best removed using a cotton wool bud wound around a round wooden stick and very little distilled water or non-aggressive solvent:

- Distilled water: Clean glass surface with the slightly moistened cotton wool bud from the center to the edge, using circular movements.
- Cleaning solution for optics, consisting of 15 % isopropanol and 85 % surgical spirit (gasoline): Clean glass surface by the slightly moistened cotton wool bud from the center to the edge, using circular movements.
- Remove dust from optical surfaces using a natural hair brush or blow using rubber bellows (air blower).
- Clean plastic parts using commercially available cleaning products (no solvents!). Tenacious impurities can be treated carefully with benzine or white spirit.
- All labels on the components may only be cleaned using a dry cotton cloth.

5.2 Maintenance

5.2.1 Wearing parts

The following wearing parts can be ordered directly from ZEISS:

Description	Order No.	Remark
Plastic plate B/W, d = 84 mm	435425-9320-000	
Glass plate, d = 84 mm	435425-9310-000	
Eyecups	444801-0000-000	required 2x
Eyepiece 10x/23		
Set of small components consisting of: <ul style="list-style-type: none"> - 2x spiral cables RJ-12 - 1x cable RJ-12 , length 1m - 2x caps, eyepiece - 1x cap, microscope body - 1x Allen key, SW 3 - 1x M6 clamping screw, microscope mount - 2x M6 clamping screws for C mount - 2x rubber rings, eyepiece - 1x Cover, C-mount Stemi 305 trino - 2x Dust cover for eyepiece tubes - 1x M4x8 screw 		
Protective dust cover, stand K	415500-1800-000	

5.3 Service

All interventions in mechanical, optical and electronic internal microscope components may only be performed by Carl Zeiss Service or by specifically **authorized** qualified staff.

To ensure your microscope is optimally set and also operates correctly over an extended period of time, we recommend that you conclude a service/maintenance contract with ZEISS.

For reordering or in case of service, contact your nearest regional ZEISS representative.

5.4 Product disposal

The product was developed, tested and produced in accordance with the valid regulations and guidelines for environmental law of the European Union.

The product and its accessories comply with the EU-regulations 2011/65/EC (RoHS) and 2012/19/EC (WEEE), and the German Electrical and Electronic Equipment Act (ElektroG).

The product contains electronic components which must be disposed of according to the requirements of the 2002/19/EC WEEE Directive and not as domestic waste. In addition, the national regulations must be complied with.

For details on disposal and recycling please refer to your relevant ZEISS sales or service organization.

6 APPENDIX**6.1 List of abbreviations**

CL	Cold Light
DIN	Deutsche Industrienorm (German Industrial Standard)
EN	Europäische Norm (European Standard)
IEC	International Electrotechnical Commission
IP	Internal Protection (Type of enclosure)
ivd	In-vitro diagnostics
LED	Light Emitting Diode
B/W	Black / white
UV	Ultra-violet

6.2 Troubleshooting

Description of fault	Cause of error	Troubleshooting
Stereo microscope illumination cannot be switched ON.	Power supply interrupted. Power switch not switched ON.	Check or establish connection to power supply. Switch power switch on back of stereo microscope ON.
	Power unit faulty.	Replace power unit acc. to operating manual Stemi 305.
Additional reflected light illuminators are not lit.	Power supply cable RJ-12 not plugged in. Reflected light not switched ON or dimmed.	Insert RJ plug correctly acc. to operating manual Stemi 305. Press and turn rotary button several times to switch reflected light ON and to increase its intensity; see operating manual Stemi 305.
Transmitted light illuminator does not light up.	Transmitted light not switched ON or dimmed.	Press and turn rotary push-button several times to switch transmitted light ON and to increase its intensity; see operating manual Stemi 305.
	Transmitted light not connected.	Open stand K EDU or K LAB and connect the plug-in connectors of the LED illumination, according to the operating manual Stemi 305.
Transmitted light illuminator dirty or control lever cannot be actuated.	Liquids or foreign matter in the transmitted light unit.	Open and stand K EDU or K LAB and eliminate foreign matter according to the operating manual Stemi 305.
Microscope mount moves down by itself.	Focusing drive too easy to move.	Set ease of motion of the focusing drive, see operating manual Stemi 305.
In mode "Connection to existing WLAN" ("Add to existing WLAN"), Stemi 305 cam cannot be connected to existing WLAN network.	Existing WLAN SSID contains a dot, e. g.: "TP-Link_02.4G".	Please rename the WLAN so that the WLAN SSID only consists of letters, digits, hyphens and underscores, e. g.: "TP-Link_02_4G".

Should any other faults occur or should you have queries regarding troubleshooting, please contact your regional ZEISS representative.

6.3 List of illustrations

Fig. 1	Warning and information labels on Stemi 305 cam	6
Fig. 2	Warning labels and apertures for LED radiation.....	6
Fig. 3	Microscope system Stemi 305 cam	8
Fig. 4	Interfaces on the Stemi 305 cam body.....	9
Fig. 5	Connecting the stereo microscope.....	16
Fig. 6	Status LED in microscope body	17
Fig. 7	Adhesive label on the back of the device	18
Fig. 8	Configuration page of Stemi 305 cam	19
Fig. 9	Safari browser of the iPad – selection of WLAN mode	20
Fig. 10	Reset of the integrated camera.....	21

6.4 Index

A

Ambient conditions 13
 Ambient temperature 13
 Apertures for LED radiation 6
 Appendix 24

C

Care 22
 Commissioning 17
 Connecting 16
 Contents 1

D

Description 8
 Device safety 3
 Dimensions 12

I

Index 27
 Indication for use 8
 Information labels 6
 Information, general 2, 16
 Installation 16
 Integrated camera 17
 Interfaces 9
 Introduction 2

L

List of abbreviations 24
 List of illustrations 26

M

Maintenance 23
 Microscope system 8

P

Product disposal 23

R

Reset 21

S

Safety 3
 Service 23
 Switching OFF 16
 Switching ON 16
 System Overview 10

T

Technical Data 12
 Troubleshooting 25

U

Use 8

W

Warning labels 6
 Warranty Notes 7
 Wearing parts 23
 Weight 12
 WLAN mode 19

Appendix to user manual

Stemi 305 cam stereo microscope

(SAP 435063-9044-008)

FCC Statement for Stemi 305 cam

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Caution:**

Changes or modifications not expressly by party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help

RF-exposure Statement for Stemi 305 cam

The device has been evaluated to meet general RF exposure requirement. The Device can be used in portable exposure conditions without restriction. Changes or modifications not expressly approved by party responsible for compliance could void the user's authority to operate the equipment.

**Caution:**

This device was tested for typical body-worn operations. To comply with R&TTE exposure requirements, a minimum separation distance of 10 mm must be maintained between the user's body, including antenna. Third-party belt-clips, holster and similar accessories used by tis device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna.

Max. SAR measurement (10 g) Head/Body: 6.2 μ W/Kg

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**Caution:**

This device was tested for typical body-worn operations. To comply with R&TTE exposure requirements, a minimum separation distance of 10 mm must be maintained between the user's body, including antenna. Third-party belt-clips, holster and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna.

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